

Anti-Cytokeratin 6 antibody (2-100) [ABT052] (STJA0039106)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | IHC/ELISA |
| Host / Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Conjugation | Unconjugated |
| Purification | The antibody was affinity-purified from ascites by affinity-chromatography using specific immunogen. |
| Dilution Range | IHC 1:200-400 ELISA 1:500-5000 |
| Formulation | PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA |
| Isotype | IgG2ak |
| Molecular Weight | Observed: 60kD |
| Storage Instruction | Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|------------------|--|
| Immunogen Region | 2-100 |
| Specificity | The antibody can specifically recognize human CK6 protein, and shows no cross reaction with CK5. |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|