

Anti-Cleaved-GSDMD antibody (N-Term) [PT0680R] (STJA0039228)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/IF/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A
Dilution Range	WB 1:2000-1:10000 IF 1:200-1:1000 ELISA 1:5000-1:20000
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Isotype	IgGk
Molecular Weight	Calculated: 53kD Observed: 35kD
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79792
Gene Symbol	GSDMD
UniProt ID	GSDMD_HUMAN
Immunogen Region	N-Term
Specificity	This antibody detects endogenous levels of Cleaved N-terminal GSDMD.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---