

Anti-Phospho-STAT5-Tyr694 antibody [PT0067R] (STJA0039276)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/IHC/IF/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A
Dilution Range	IHC 1:100-1:500 WB 1:2000-1:10000 IF 1:200-1:1000 ELISA 1:5000-1:20000
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Isotype	IgGk
Molecular Weight	Calculated: 91kD Observed: 105kD
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Specificity	This antibody detects endogenous levels of STAT5A/B only when phosphorylated at Tyr694.This antibody does not recognize phosphorylated at other sites.
-------------	--

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---