

**Anti-Pan-Akt antibody [PT0654R] (STJA0039389)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/IHC/IF/IP/ELISA |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse/Rat    |

**PRODUCT PROPERTIES**

|                     |   |
|---------------------|---|
| Clonality           | Monoclonal  |
| Conjugation         | Unconjugated  |
| Purification        | Protein A   |
| Dilution Range      | IHC 1:200-1:1000<br>WB 1:2000-1:10000<br>IF 1:200-1:1000<br>ELISA 1:5000-1:20000<br>IP 1:50-1:200 |
| Formulation         | PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA  |
| Isotype             | IgGk  |
| Molecular Weight    | Calculated: 55kD<br>Observed: 55kD  |
| Storage Instruction | Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.     |

**TARGET INFORMATION**

|             |            |
|-------------|------------|
| Specificity | Endogenous |
|-------------|------------|

**ADDITIONAL INFORMATION**

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|