

Anti-Cleaved-Caspase-3-Asp175 antibody [PT0457R] (STJA0039448)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/IHC/IF/IP/ELISA
	Host / Source	Rabbit
	Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Conjugation	Unconjugated
	Purification	Protein A
	Dilution Range	IHC 1:200-1:1000 WB 1:1000-1:5000 IF 1:200-1:1000 ELISA 1:5000-1:20000 IP 1:50-1:200
	Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
	Isotype	IgGk
	Molecular Weight	Calculated: 17kD, 19kD Observed: 17kD, 19kD
	Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Specificity	This antibody detects endogenous levels of Cleaved Caspase-3 (Asp175).
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.