

**Anti-NF-Kappa B1 p105/p50 antibody [PT0463R] (STJA0039451)**
**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/IHC/IF/IP/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A
Dilution Range	IHC 1:200-1:1000 WB 1:2000-1:10000 IF 1:200-1:1000 ELISA 1:5000-1:20000 IP 1:50-1:200
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Isotype	IgGk
Molecular Weight	Calculated: 50kD, 105kD Observed: 50kD, 120kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

**TARGET INFORMATION**

Specificity	Endogenous
-------------	------------

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---