

**Anti-Histone H3 antibody [PT0508R] (STJA0039467)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/IHC/IF/IP/ELISA/CHIP/Cut&Tag
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A
Dilution Range	IHC 1:5000-1:20000 WB 1:5000-1:20000 IF 1:200-1:1000 ELISA 1:5000-1:20000 IP 1:50-1:200 CHIP 1:50-1:100 Cut&Tag 1:50-1:100
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Isotype	IgGk
Molecular Weight	Calculated: 15kD Observed: 15kD
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

**TARGET INFORMATION**

Specificity	Endogenous
-------------	------------

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---