

Anti-Phospho-Rpb1-Ser5 antibody [PT0923R] (STJA0039821)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	WB/IHC/IF/IP/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES	
Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A
Dilution Range	IHC 1:200-1:1000 WB 1:500-1:2000 IF 1:200-1:1000 ELISA 1:5000-1:20000 IP 1:50-1:200
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Isotype	IgGk
Molecular Weight	Calculated: 217kD Observed: 217kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	5430
Gene Symbol	POLR2A
UniProt ID	RPB1_HUMAN
Specificity	Rpb1 (Phospho Ser5) Monoclonal Antibody detects endogenous levels of Rpb1 protein only when phosphorylated at Ser5. The name of modified sites may be influenced by many factors, such as species (the modified site was not originally found in human samp
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.