

## Anti-Phospho-TrkA/TrkB-Tyr516 antibody [PT1381R] (STJA0040021)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/IP/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A
Dilution Range	WB 1:2000-1:10000 IF 1:200-1:1000 ELISA 1:5000-1:20000 IP 1:50-1:200
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Isotype	IgGk
Molecular Weight	Calculated: 92kD, 87kD Observed: 92kD, 87kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Specificity	TrkA (Phospho Tyr490) /TrkB (Phospho Tyr516) Antibody detects endogenous levels of TrkA/TrkB protein only when phosphorylated at Y490/Y516. The name of modified sites may be influenced by many factors, such as species (the modified site was not origin
-------------	---

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---