

## Anti-ERCC1 antibody [PT1303R] (STJA0040098)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/IHC/IF/ELISA    |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse/Rat    |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Conjugation         | Unconjugated   |
| Purification        | Protein A  |
| Dilution Range      | IHC 1:200-1:1000<br>WB 1:2000-1:10000<br>IF 1:200-1:1000<br>ELISA 1:5000-1:20000               |
| Formulation         | PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA   |
| Isotype             | IgGk   |
| Molecular Weight    | Calculated: 33kD<br>Observed: 39kD   |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">2067</a>        |
| Gene Symbol | <a href="#">ERCC1</a>       |
| UniProt ID  | <a href="#">ERCC1_HUMAN</a> |
| Specificity | Endogenous                  |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|