

Anti-Phospho-HER2/HER4-Tyr1284 antibody [PT1236R] (STJA0040156)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/IF/IP/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A
Dilution Range	WB 1:2000-1:10000 IF 1:200-1:1000 ELISA 1:5000-1:20000 IP 1:50-1:200
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Isotype	IgGk
Molecular Weight	Calculated: 147kD Observed: 180kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Specificity	HER2 (Phospho Tyr1248) /HER4 (Phospho Tyr1284) Antibody detects endogenous levels of HER2/HER4 protein only when phosphorylated at Y1248 (HER2) /Tyr1284 (HER4).The name of modified sites may be influenced by many factors, such as species (the modifie
-------------	--

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---