

Anti-Phospho-c-Myb-Ser11 antibody [PT1200R] (STJA0040193)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/IHC/IF/IP/ELISA
	Host / Source	Rabbit
	Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Conjugation	Unconjugated
	Purification	Protein A
	Dilution Range	IHC 1:200-1:1000 WB 1:2000-1:10000 IF 1:200-1:1000 ELISA 1:5000-1:20000 IP 1:50-1:200
	Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
	Isotype	IgGk
	Molecular Weight	Calculated: 72kD Observed: 80kD
	Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Specificity	c-Myb (Phospho Ser11) Antibody detects endogenous levels of c-Myb protein only when phosphorylated at S11.The name of modified sites may be influenced by many factors, such as species (the modified site was not originally found in human samples) and
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.