

Anti-CD127 antibody [PTF0044] (STJA0040310)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IF/WB
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A
Dilution Range	WB 1:4000-1:20000 IF 1:200-1:1000 FC: 1 µg per million cells in 100 µl volume. For optimal results, the reagent should be titrated according to each specific application
Formulation	PBS, pH 7.4
Isotype	IgG1k
Molecular Weight	Calculated: 52kD Observed: 70kD
Storage Instruction	Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

TARGET INFORMATION

Specificity	Endogenous
-------------	------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---