

Anti-CD5 antibody [PTF0040] (STJA0040313)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IF/WB
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A
Dilution Range	WB 1:1000-1:5000 IF 1:200-1:1000 FC: 1 µg per million cells in 100 µl volume. For optimal results, the reagent should be titrated according to each specific application
Formulation	PBS, pH 7.4
Isotype	IgG1k
Molecular Weight	Calculated: 55kD Observed: 67kD
Storage Instruction	Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

TARGET INFORMATION

Gene ID	921
Gene Symbol	CD5
UniProt ID	CD5_HUMAN
Specificity	Endogenous

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---