

Anti-CD79b antibody [PTF0050] (STJA0040314)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IF/WB/IHC
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A
Dilution Range	IHC 1:100-1:200 WB 1:2000-1:10000 IF 1:200-1:1000 FC:1 µg per million cells in 100 µl volume. For optimal results, the reagent should be titrated according to each specific application
Formulation	PBS, pH 7.4
Isotype	IgG1k
Molecular Weight	Calculated: 26kD Observed: 24kD
Storage Instruction	Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.

TARGET INFORMATION

Gene ID	974
Gene Symbol	CD79B
UniProt ID	CD79B_HUMAN
Specificity	Endogenous

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---