

## Anti-Phospho-ABL1-Tyr245 antibody [PT1402R] (STJA0040363)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/ELISA
Host / Source	Rabbit
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A
Dilution Range	WB 1:2000-1:10000 IF 1:200-1:1000 ELISA 1:5000-1:20000
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Isotype	IgGk
Molecular Weight	Calculated: 123kD Observed: 210kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">25</a>
Gene Symbol	<a href="#">ABL1</a>
UniProt ID	<a href="#">ABL1_HUMAN</a>
Specificity	ABL1 (Phospho Tyr245) Antibody detects endogenous levels of ABL1 protein only when phosphorylated at Y245. The name of modified sites may be influenced by many factors, such as species (the modified site was not originally found in human samples) and

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---