

Anti-Tyrosine Hydroxylase antibody [LNC1] (STJA0040396)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/ICC
Host / Source	Mouse
Reactivity	All Mammalian

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/ml
Conjugation	Unconjugated
Purification	Protein G purified culture supernatant.
Dilution Range	WB 1:1000 IHC 1:1000 ICC 1:500-1:1000
Formulation	100 µl in 10 mM Tris, 50 mM NaCl, 0.065% Sodium Azide pH 7.4
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	25085
Gene Symbol	Th
UniProt ID	TY3H_RAT
Specificity	Specific for endogenous levels of the ~60 kDa tyrosine hydroxylase protein.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---