

**Anti-Recombinant-PSD 95 antibody (77-299 aa) [K28/43] {FL490 Ex: 491 nm, Em: 515 nm}
(STJA0040405)**

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC
Host / Source	Chicken
Reactivity	Chicken/Human/Mouse/Non-Human Primates/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	FL490 Ex: 491 nm, Em: 515 nm
Purification	Purified by affinity chromatography
Dilution Range	IHC 1:500 ICC 1:500
Formulation	1X PBS, 0.05% Sodium Azide pH 7.4
Isotype	IgY
Molecular Weight	95-110 kDa (varies with cell background due to phosphorylation)
Storage Instruction	Aliquot and store at ≤-20°C for long term storage. For short term storage, store at 2-8°C. For maximum recovery of product, centrifuge the vial prior to removing the cap.

TARGET INFORMATION

Gene ID	1742
Gene Symbol	DLG4
UniProt ID	DLG4_HUMAN
Immunogen Region	77-299 aa
Immunogen Sequence	Human
Specificity	No cross-reactivity against recombinant Chapsyn-110, SAP97 or SAP102

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---