

Anti-Recombinant-VGlut1 antibody (C-Term) [N28/9] (STJA0040417)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/IHC/ICC/FC/ELISA
Host / Source	Chicken
Reactivity	Fish/Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by affinity chromatography
Dilution Range	WB 1:500-1000 IHC 1:100-1:200 ICC 1:100-1:200
Formulation	PBS
Isotype	IgY
Molecular Weight	52 kDa
Storage Instruction	Aliquot and store at $\leq -20^{\circ}\text{C}$ for long term storage. For short term storage, store at $2-8^{\circ}\text{C}$. For maximum recovery of product, centrifuge the vial prior to removing the cap.

TARGET INFORMATION

Gene ID	116638
Gene Symbol	Slc17a7
UniProt ID	VGLU1_RAT
Immunogen Region	C-Term
Immunogen Sequence	Rat
Specificity	Does not cross-react with VGlut2 or VGlut3

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---