

**Anti-Recombinant-RGS14 antibody (Full-Length) [N133/21] (STJA0040442)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/IHC/ICC/ELISA
Host / Source	Chicken
Reactivity	Mouse/Rat

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by affinity chromatography
Formulation	PBS
Isotype	IgY
Molecular Weight	60 kDa
Storage Instruction	Aliquot and store at $\leq 20^{\circ}\text{C}$ for long term storage. For short term storage, store at $2-8^{\circ}\text{C}$ . For maximum recovery of product, centrifuge the vial prior to removing the cap.

**TARGET INFORMATION**

Gene ID	<a href="#">114705</a>
Gene Symbol	<a href="#">Rgs14</a>
UniProt ID	<a href="#">RGS14_RAT</a>
Immunogen Region	Full-Length
Immunogen Sequence	Rat
Specificity	Does not react with RGS12

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------