

**Anti-Recombinant-RGS14 antibody (Full-Length) [N133/21] {FL490 Ex: 491 nm, Em: 515 nm}
(STJA0040443)**

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/ICC
Host / Source	Chicken
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	FL490 Ex: 491 nm, Em: 515 nm
Purification	Purified by affinity chromatography
Formulation	1X PBS, 0.05% Sodium Azide pH 7.4
Isotype	IgY
Molecular Weight	60 kDa
Storage Instruction	Aliquot and store at $\leq -20^{\circ}\text{C}$ for long term storage. For short term storage, store at $2-8^{\circ}\text{C}$. For maximum recovery of product, centrifuge the vial prior to removing the cap.

TARGET INFORMATION

Gene ID	114705
Gene Symbol	Rgs14
UniProt ID	RGS14_RAT
Immunogen Region	Full-Length
Immunogen Sequence	Rat
Specificity	Does not react with RGS12

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---