

Anti-Recombinant-Olig1 antibody (N-Term) [N149/25] (STJA0040448)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/ICC/ELISA
Host / Source	Chicken
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by affinity chromatography
Formulation	PBS
Isotype	IgY
Molecular Weight	30 kDa
Storage Instruction	Aliquot and store at $\leq 20^{\circ}\text{C}$ for long term storage. For short term storage, store at $2-8^{\circ}\text{C}$. For maximum recovery of product, centrifuge the vial prior to removing the cap.

TARGET INFORMATION

Gene ID	50914
Gene Symbol	Olig1
UniProt ID	OLIG1_MOUSE
Immunogen Region	N-Term
Immunogen Sequence	Mouse
Specificity	No cross-reactivity against Olig2

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---