

Anti-Recombinant-Dardarin/LRRK2, C-terminus antibody (C-Term) [N241A/34] (STJA0040460)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/ICC/IP/ELISA
Host / Source	Chicken
Reactivity	Drosophila/Human/Mouse/Non-Human Primates/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by affinity chromatography
Dilution Range	WB 1:1000
Formulation	1X PBS, 0.05% Sodium Azide pH 7.4
Isotype	IgY
Molecular Weight	>200 kDa
Storage Instruction	Aliquot and store at $\leq -20^{\circ}\text{C}$ for long term storage. For short term storage, store at $2-8^{\circ}\text{C}$. For maximum recovery of product, centrifuge the vial prior to removing the cap.

TARGET INFORMATION


Gene ID	120892
Gene Symbol	LRRK2
UniProt ID	LRRK2_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	Human
Specificity	No cross-reactivity reported

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



IHC staining rat brain cerebral cortex using anti-LRRK2 antibody green fluorescence STJA0040460



IHC staining rat brain cerebral cortex with anti-LRRK2 antibody green and DAPI blue STJA0040460



Western blot anti-LRRK2 chicken Mab detection in LRRK2-FLAG transfected COS-7 lysates STJA0040460