

**Anti-Recombinant-Dardarin/LRRK2, C-terminus antibody (C-Term) [N241A/34] (STJA0040461)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/IHC/ICC/IP/ELISA
Host / Source	Chicken
Reactivity	Drosophila/Human/Mouse/Non-Human Primates/Rat

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by affinity chromatography
Dilution Range	WB 1:1000
Formulation	PBS
Isotype	IgY
Molecular Weight	>200 kDa
Storage Instruction	Aliquot and store at $\leq -20^{\circ}\text{C}$ for long term storage. For short term storage, store at $2-8^{\circ}\text{C}$ . For maximum recovery of product, centrifuge the vial prior to removing the cap.

**TARGET INFORMATION**

Gene ID	<a href="#">120892</a>
Gene Symbol	<a href="#">LRRK2</a>
UniProt ID	<a href="#">LRRK2_HUMAN</a>
Immunogen Region	C-Term
Immunogen Sequence	Human
Specificity	No cross-reactivity reported

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---