

Anti-Recombinant-Arl13b antibody (C-Term) [N295B/66] (STJA0040467)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/ICC/ELISA
Host / Source	Chicken
Reactivity	Human/Mouse/Non-Human Primates/Rat/Zebrafish

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by affinity chromatography
Dilution Range	WB 1:1000
Formulation	1X PBS, 0.05% Sodium Azide pH 7.4
Isotype	IgY
Molecular Weight	50 kDa
Storage Instruction	Aliquot and store at $\leq 20^{\circ}\text{C}$ for long term storage. For short term storage, store at $2-8^{\circ}\text{C}$. For maximum recovery of product, centrifuge the vial prior to removing the cap.

TARGET INFORMATION


Gene ID	68146
Gene Symbol	Arl13b
UniProt ID	AR13B_MOUSE
Immunogen Region	C-Term
Immunogen Sequence	Mouse
Specificity	Does not cross-react with Arl13a (based on KO validation results)

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



IHC staining rat brain lateral ventricle with anti-Arl13b antibody ciliated cells in green STJA0040467



IHC staining rat brain lateral ventricle with anti-Arl13b antibody ciliated cells STJA0040467