

Anti-Recombinant-Beta4-spectrin antibody (C-Term) [N393/76] {FL490 Ex: 491 nm, Em: 515 nm} (STJA0040476)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC
Host / Source	Chicken
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	FL490 Ex: 491 nm, Em: 515 nm
Purification	Purified by affinity chromatography
Formulation	1X PBS, 0.05% Sodium Azide pH 7.4
Isotype	IgY
Molecular Weight	>200 kDa
Storage Instruction	Aliquot and store at $\leq -20^{\circ}\text{C}$ for long term storage. For short term storage, store at $2-8^{\circ}\text{C}$. For maximum recovery of product, centrifuge the vial prior to removing the cap.

TARGET INFORMATION

Gene ID	57731
Gene Symbol	SPTBN4
UniProt ID	SPTN4_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	Human
Specificity	Does not cross-react with other Beta-spectrins

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---