

Anti-Recombinant-VACHT antibody (55-125 aa) [N425/45] (STJA0040481)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB/IHC/ICC/ELISA
Host / Source	Chicken
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by affinity chromatography
Formulation	PBS
Isotype	IgY
Molecular Weight	60 kDa
Storage Instruction	Aliquot and store at $\leq 20^{\circ}\text{C}$ for long term storage. For short term storage, store at $2-8^{\circ}\text{C}$. For maximum recovery of product, centrifuge the vial prior to removing the cap.

TARGET INFORMATION

Gene ID	60422
Gene Symbol	Slc18a3
UniProt ID	VACHT_RAT
Immunogen Region	55-125 aa
Immunogen Sequence	Rat
Specificity	Does not cross-react with VMAT1 or VMAT2

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---