

Anti-Recombinant-Homer1 antibody (C-Term) [L113/130] {FL550 Ex: 550 nm, Em: 575 nm} (STJA0040495)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/ICC/AT
Host / Source	Chicken
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	FL550 Ex: 550 nm, Em: 575 nm
Purification	Purified by affinity chromatography
Dilution Range	IHC 1:100 ICC 1:100
Formulation	1X PBS, 0.05% Sodium Azide pH 7.4
Isotype	IgY
Molecular Weight	50 kDa
Storage Instruction	Aliquot and store at $\leq -20^{\circ}\text{C}$ for long term storage. For short term storage, store at $2-8^{\circ}\text{C}$. For maximum recovery of product, centrifuge the vial prior to removing the cap.

TARGET INFORMATION

Gene ID	26556
Gene Symbol	Homer1
UniProt ID	HOME1_MOUSE
Immunogen Region	C-Term
Immunogen Sequence	Mouse
Specificity	No cross-reactivity reported

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---