

Anti-Recombinant-VGAT antibody (N-Term) [L118/80] {FL490 Ex: 491 nm, Em: 515 nm} (STJA0040500)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC
Host / Source	Chicken
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	FL490 Ex: 491 nm, Em: 515 nm
Purification	Purified by affinity chromatography
Dilution Range	IHC 1:100 ICC 1:100
Formulation	1X PBS, 0.05% Sodium Azide pH 7.4
Isotype	IgY
Molecular Weight	50 kDa
Storage Instruction	Aliquot and store at $\leq -20^{\circ}\text{C}$ for long term storage. For short term storage, store at 2-8°C. For maximum recovery of product, centrifuge the vial prior to removing the cap.

TARGET INFORMATION

Gene ID	22348
Gene Symbol	Slc32a1
UniProt ID	VIAAT_MOUSE
Immunogen Region	N-Term
Immunogen Sequence	Mouse
Specificity	No cross-reactivity reported

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---