

## Anti-Recombinant-Parvalbumin antibody (Full-Length) [L114/81] {FL490 Ex: 491 nm, Em: 515 nm} (STJA0040506)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC
Host / Source	Chicken
Reactivity	Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	FL490 Ex: 491 nm, Em: 515 nm
Purification	Purified by affinity chromatography
Formulation	1X PBS, 0.05% Sodium Azide pH 7.4
Isotype	IgY
Molecular Weight	12 kDa
Storage Instruction	Aliquot and store at $\leq -20^{\circ}\text{C}$ for long term storage. For short term storage, store at $2-8^{\circ}\text{C}$ . For maximum recovery of product, centrifuge the vial prior to removing the cap.

### TARGET INFORMATION

Immunogen Region	Full-Length
Immunogen Sequence	Rat
Specificity	No cross-reactivity reported

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---