

**Anti-Recombinant-Calretinin antibody (1-178 N-Term) [L122/6] (STJA0040511)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | ELISA/ICC/IHC/WB   |
| Host / Source | Chicken            |
| Reactivity    | Human/Mouse/Rat    |

**PRODUCT PROPERTIES**

|                     |                                                                                                                                                                                                                |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                                                                                                                                     |
| Concentration       | 1 mg/mL                                                                                                                                                                                                        |
| Conjugation         | Unconjugated                                                                                                                                                                                                   |
| Purification        | Purified by affinity chromatography                                                                                                                                                                            |
| Formulation         | PBS                                                                                                                                                                                                            |
| Isotype             | IgY                                                                                                                                                                                                            |
| Molecular Weight    | 25-35 kDa                                                                                                                                                                                                      |
| Storage Instruction | Aliquot and store at $\leq 20^{\circ}\text{C}$ for long term storage. For short term storage, store at $2-8^{\circ}\text{C}$ . For maximum recovery of product, centrifuge the vial prior to removing the cap. |

**TARGET INFORMATION**

|                    |                              |
|--------------------|------------------------------|
| Gene ID            | <a href="#">794</a>          |
| Gene Symbol        | <a href="#">CALB2</a>        |
| UniProt ID         | <a href="#">CALB2_HUMAN</a>  |
| Immunogen Region   | 1-178 N-Term                 |
| Immunogen Sequence | Human                        |
| Specificity        | No cross-reactivity reported |

**ADDITIONAL INFORMATION**

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|