

Anti-Recombinant-Bassoon antibody (157-517 aa) [L124/59] {FL490 Ex: 491 nm, Em: 515 nm} (STJA0040518)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	IHC/ICC/AT
Host / Source	Chicken
Reactivity	Mouse/Rat
PRODUCT PROPERTIES	
Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	FL490 Ex: 491 nm, Em: 515 nm
Purification	Purified by affinity chromatography
Dilution Range	IHC 1:100 ICC 1:100
Formulation	1X PBS, 0.05% Sodium Azide pH 7.4
Isotype	IgY
Molecular Weight	400 kDa for full-length protein, many smaller bands possible depending on sample used
Storage Instruction	Aliquot and store at ≤-20°C for long term storage. For short term storage, store at 2-8°C. For maximum recovery of product, centrifuge the vial prior to removing the cap.
TARGET INFORMATION	
Gene Symbol	Bsn Uniprot ID= BSN_RAT "
Immunogen Region	157-517 aa
Immunogen Sequence	Rat
Specificity	No cross-reactivity reported
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.