

Anti-Recombinant-c-FOS antibody (Full-Length) [N486/32] (STJA0040522)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/ICC/ELISA
Host / Source	Chicken
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by affinity chromatography
Formulation	1X PBS, 0.05% Sodium Azide pH 7.4
Isotype	IgY
Molecular Weight	70 kDa
Storage Instruction	Aliquot and store at $\leq -20^{\circ}\text{C}$ for long term storage. For short term storage, store at $2-8^{\circ}\text{C}$. For maximum recovery of product, centrifuge the vial prior to removing the cap.

TARGET INFORMATION


Gene ID	2353
Gene Symbol	FOS
UniProt ID	FOS_HUMAN
Immunogen Region	Full-Length
Immunogen Sequence	Human
Specificity	No cross-reactivity reported

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




anti-c-FOS antibody IHC staining in rat brain sections neuronal nuclei STJA0040522




Western blot anti-c-Fos antibody detection c-Fos protein in HeLa cells at 65 kDa STJA0040522



ICC staining HeLa cells with anti-c-Fos antibody nuclear localization in green STJA0040522



IHC staining rat hippocampus with anti-c-FOS antibody neuronal nuclei STJA0040522



IHC staining rat hippocampus using anti-c-Fos antibody in seizure vs control tissue STJA0040522