

Anti-Recombinant-GAD65/67 antibody (Full-Length) [L127/12] (STJA0040529)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | WB/IHC/ICC/ELISA |
| Host / Source | Chicken |
| Reactivity | Human/Mouse/Rat |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Concentration | 1 mg/mL |
| Conjugation | Unconjugated |
| Purification | Purified by affinity chromatography |
| Dilution Range | WB 1:1000 IHC 1:100 ICC 1:100 |
| Formulation | PBS |
| Isotype | IgY |
| Molecular Weight | 70 kDa |
| Storage Instruction | Aliquot and store at $\leq -20^{\circ}\text{C}$ for long term storage. For short term storage, store at 2-8°C. For maximum recovery of product, centrifuge the vial prior to removing the cap. |

TARGET INFORMATION

| | |
|--------------------|---|
| Gene ID | 2571 |
| Gene Symbol | GAD1 |
| UniProt ID | DCE1_HUMAN |
| Immunogen Region | Full-Length |
| Immunogen Sequence | Human |
| Specificity | Cross-reacts with GAD25, GAD44 and GAD65/GAD2 |

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|