

**Anti-Recombinant-GAD65/67 antibody (Full-Length) [L127/12] {FL550 Ex: 550 nm, Em: 575 nm}
(STJA0040531)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ICC/IHC
Host / Source	Chicken
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	FL550 Ex: 550 nm, Em: 575 nm
Purification	Purified by affinity chromatography
Dilution Range	IHC 1:100 ICC 1:100
Formulation	1X PBS, 0.05% Sodium Azide pH 7.4
Isotype	IgY
Molecular Weight	70 kDa
Storage Instruction	Aliquot and store at $\leq -20^{\circ}\text{C}$ for long term storage. For short term storage, store at 2-8°C. For maximum recovery of product, centrifuge the vial prior to removing the cap.

TARGET INFORMATION

Gene ID	2571
Gene Symbol	GAD1
UniProt ID	DCE1_HUMAN
Immunogen Region	Full-Length
Immunogen Sequence	Human
Specificity	Cross-reacts with GAD25, GAD44 and GAD65/GAD2

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---