

**Anti-Recombinant-ALB, Human Serum Albumin antibody [3B8-H3-G1] {FL490 Ex: 491 nm, Em: 515 nm} (STJA0040542)**

## GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC/ICC            |
| Host / Source | Chicken            |
| Reactivity    | Human              |

## PRODUCT PROPERTIES

|                     |                                                                                                                                                                                                                 |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                                                                                                                                      |
| Concentration       | 0.5 mg/mL                                                                                                                                                                                                       |
| Conjugation         | FL490 Ex: 491 nm, Em: 515 nm                                                                                                                                                                                    |
| Purification        | Purified by affinity chromatography                                                                                                                                                                             |
| Dilution Range      | ICC 1:100                                                                                                                                                                                                       |
| Formulation         | 1X PBS, 0.05% Sodium Azide pH 7.4                                                                                                                                                                               |
| Isotype             | IgY                                                                                                                                                                                                             |
| Molecular Weight    | 66 kDa                                                                                                                                                                                                          |
| Storage Instruction | Aliquot and store at $\leq -20^{\circ}\text{C}$ for long term storage. For short term storage, store at $2-8^{\circ}\text{C}$ . For maximum recovery of product, centrifuge the vial prior to removing the cap. |

## TARGET INFORMATION

|                    |                              |
|--------------------|------------------------------|
| Gene ID            | <a href="#">213</a>          |
| Gene Symbol        | <a href="#">ALB</a>          |
| UniProt ID         | <a href="#">ALBU_HUMAN</a>   |
| Immunogen Sequence | Human                        |
| Specificity        | No cross-reactivity reported |

## ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|