

Anti-Recombinant-ALB, Human Serum Albumin antibody [3B8-H3-G1] {FL550 Ex: 550 nm, Em: 575 nm} (STJA0040543)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/ICC
Host / Source	Chicken
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	FL550 Ex: 550 nm, Em: 575 nm
Purification	Purified by affinity chromatography
Dilution Range	ICC 1:100
Formulation	1X PBS, 0.05% Sodium Azide pH 7.4
Isotype	IgY
Molecular Weight	66 kDa
Storage Instruction	Aliquot and store at $\leq -20^{\circ}\text{C}$ for long term storage. For short term storage, store at $2-8^{\circ}\text{C}$. For maximum recovery of product, centrifuge the vial prior to removing the cap.

TARGET INFORMATION

Gene ID	213
Gene Symbol	ALB
UniProt ID	ALBU_HUMAN
Immunogen Sequence	Human
Specificity	No cross-reactivity reported

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---