

## Anti-Recombinant-Cortactin antibody [4F11-G5-F10] (STJA0040546)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/ICC/IP
Host / Source	Chicken
Reactivity	Bovine/Chicken/Hamster/Human/Mouse/Rat

### PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by affinity chromatography
Dilution Range	WB 1:1000-1:3000 IHC 1:100-1:500 ICC 1:100-1:500
Formulation	1X PBS, 0.05% Sodium Azide pH 7.4
Isotype	IgY
Molecular Weight	65 kDa
Storage Instruction	Aliquot and store at $\leq -20^{\circ}\text{C}$ for long term storage. For short term storage, store at 2-8 $^{\circ}\text{C}$ . For maximum recovery of product, centrifuge the vial prior to removing the cap.

### TARGET INFORMATION


Immunogen Sequence	Chicken
Specificity	No cross-reactivity reported

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



ICC staining HeLa cells with anti-Cortactin antibody green fluorescence and DAPI nuclear stain STJA0040546



Western blot WB HeLa cell lysates using anti-Cortactin recombinant chicken mAb at 65 kDa STJA0040546