

**Anti-Recombinant-Cortactin antibody [4F11-G5-F10] {FL490 Ex: 491 nm, Em: 515 nm}
(STJA0040548)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IHC/ICC
Host / Source	Chicken
Reactivity	Bovine/Chicken/Hamster/Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	FL490 Ex: 491 nm, Em: 515 nm
Purification	Purified by affinity chromatography
Dilution Range	ICC 1:100-1:500
Formulation	1X PBS, 0.05% Sodium Azide pH 7.4
Isotype	IgY
Molecular Weight	65 kDa
Storage Instruction	Aliquot and store at $\leq -20^{\circ}\text{C}$ for long term storage. For short term storage, store at $2-8^{\circ}\text{C}$. For maximum recovery of product, centrifuge the vial prior to removing the cap.

TARGET INFORMATION

Immunogen Sequence	Chicken
Specificity	No cross-reactivity reported

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---