

Anti-Recombinant-Cytokeratin 8 antibody (472-483) [E8-C6-B4] (STJA0040552)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/ICC
Host / Source	Chicken
Reactivity	Bovine/Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by affinity chromatography
Dilution Range	WB 1:1000-1:3000 IHC 1:500-1:2000 ICC 1:500-1:2000
Formulation	1X PBS, 0.05% Sodium Azide pH 7.4
Isotype	IgY
Molecular Weight	52 kDa
Storage Instruction	Aliquot and store at $\leq -20^{\circ}\text{C}$ for long term storage. For short term storage, store at $2-8^{\circ}\text{C}$. For maximum recovery of product, centrifuge the vial prior to removing the cap.

TARGET INFORMATION


Gene ID	3856
Gene Symbol	KRT8
UniProt ID	K2C8_HUMAN
Immunogen Region	472-483
Immunogen Sequence	Human
Specificity	No cross-reactivity reported

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




ICC staining HeLa cells with anti-Cytokeratin 8 antibody cytoskeleton in green STJA0040552



Western blot anti-Cytokeratin 8 antibody detection at 52 kDa in HeLa cell lysates STJA0040552



ICC staining HeLa cells with anti-Cytokeratin 8 antibody cytoskeleton in green STJA0040552



ICC staining HeLa cells using anti-Cytokeratin 8 antibody green cytoskeleton with DAPI STJA0040552