

Anti-Recombinant-Cytokeratin 8 antibody (472-483) [E8-C6-B4] {FL650 Ex: 655 nm, Em: 676 nm} (STJA0040557)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/ICC
Host / Source	Chicken
Reactivity	Bovine/Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	FL650 Ex: 655 nm, Em: 676 nm
Purification	Purified by affinity chromatography
Dilution Range	ICC 1:500-1:2000
Formulation	1X PBS, 0.05% Sodium Azide pH 7.4
Isotype	IgY
Molecular Weight	52 kDa
Storage Instruction	Aliquot and store at $\leq -20^{\circ}\text{C}$ for long term storage. For short term storage, store at $2-8^{\circ}\text{C}$. For maximum recovery of product, centrifuge the vial prior to removing the cap.

TARGET INFORMATION

Gene ID	3856
Gene Symbol	KRT8
UniProt ID	K2C8_HUMAN
Immunogen Region	472-483
Immunogen Sequence	Human
Specificity	No cross-reactivity reported

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---