

## Anti-Recombinant-M6PR antibody [2G11] (STJA0040559)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/ICC/IP
Host / Source	Chicken
Reactivity	Bovine/Human/Non-Human Primates

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by affinity chromatography
Dilution Range	IHC 1:100-1:500 ICC 1:100-1:500
Formulation	PBS
Isotype	IgY
Molecular Weight	300 kDa
Storage Instruction	Aliquot and store at $\leq -20^{\circ}\text{C}$ for long term storage. For short term storage, store at $2-8^{\circ}\text{C}$ . For maximum recovery of product, centrifuge the vial prior to removing the cap.

### TARGET INFORMATION

Immunogen Sequence	Bovine
Specificity	No cross-reactivity reported

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---