

Anti-Recombinant-M6PR antibody [2G11] {FL550 Ex: 550 nm, Em: 575 nm} (STJA0040561)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	IHC/ICC
Host / Source	Chicken
Reactivity	Bovine/Human/Non-Human Primates
PRODUCT PROPERTIES	
Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	FL550 Ex: 550 nm, Em: 575 nm
Purification	Purified by affinity chromatography
Dilution Range	IHC 1:100-1:500 ICC 1:100-1:500
Formulation	1X PBS, 0.05% Sodium Azide pH 7.4
Isotype	IgY
Molecular Weight	300 kDa
Storage Instruction	Aliquot and store at ≤-20°C for long term storage. For short term storage, store at 2-8°C. For maximum recovery of product, centrifuge the vial prior to removing the cap.
TARGET INFORMATION	
Immunogen Sequence	Bovine
Specificity	No cross-reactivity reported
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.