

Anti-Recombinant-Rab 9 antibody [Clone 9] (STJA0040565)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	WB/IHC/ICC/IP
Host / Source	Chicken
Reactivity	Bovine/Canine/Hamster/Human/Mouse/Non-Human Primates/Rat
PRODUCT PROPERTIES	
Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified by affinity chromatography
Dilution Range	WB 1:1000 IHC 1:50-1:200 ICC 1:50-1:200
Formulation	PBS
Isotype	IgY
Molecular Weight	23 kDa
Storage Instruction	Aliquot and store at $\leq -20^{\circ}\text{C}$ for long term storage. For short term storage, store at $2-8^{\circ}\text{C}$. For maximum recovery of product, centrifuge the vial prior to removing the cap.
TARGET INFORMATION	
Immunogen Sequence	Canine
Specificity	No cross-reactivity reported
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.