

**Anti-Recombinant-Rab 9 antibody [Clone 9] {FL490 Ex: 491 nm, Em: 515 nm} (STJA0040566)**
**GENERAL INFORMATION**

|               |  |
|---------------|--|
| Product Type  | Primary antibodies                                       |
| Applications  | IHC/ICC  |
| Host / Source | Chicken  |
| Reactivity    | Bovine/Canine/Hamster/Human/Mouse/Non-Human Primates/Rat |

**PRODUCT PROPERTIES**

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Concentration       | 0.5 mg/mL  |
| Conjugation         | FL490 Ex: 491 nm, Em: 515 nm   |
| Purification        | Purified by affinity chromatography  |
| Dilution Range      | IHC 1:50-1:200<br>ICC 1:50-1:200   |
| Formulation         | 1X PBS, 0.05% Sodium Azide pH 7.4  |
| Isotype             | IgY  |
| Molecular Weight    | 23 kDa   |
| Storage Instruction | Aliquot and store at $\leq -20^{\circ}\text{C}$ for long term storage. For short term storage, store at 2-8°C. For maximum recovery of product, centrifuge the vial prior to removing the cap. |

**TARGET INFORMATION**

|                    |                              |
|--------------------|------------------------------|
| Immunogen Sequence | Canine                       |
| Specificity        | No cross-reactivity reported |

**ADDITIONAL INFORMATION**

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|