

Anti-Actin antibody (STJA0040577)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Chicken
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	10 mg/mL
Conjugation	Unconjugated
Purification	Mixture of IgY fraction and affinity-purified antibodies
Dilution Range	WB 1:1, 000-1:2, 000 WB Brain 1:1, 000-1:2, 000
Formulation	Sodium phosphate (10 mM, pH 7.2) buffered isotonic saline (0.9%, w/v) , glycerol (50%, v/v) , with sodium azide (0.02%, w/v) as an anti-microbial agent.
Isotype	IgY
Molecular Weight	42 kDa
Storage Instruction	Store at -20°C in the dark. Under these conditions, the antibodies should have a shelf life of at least twelve months, provided they remain sterile.

TARGET INFORMATION

Gene ID	72
Gene Symbol	ACTG2
UniProt ID	ACTH_HUMAN
Immunogen Sequence	Rabbit

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot using Actin as loading control for TFPI protein expression in HeLa cells STJA0040577