

Anti-Beta-Gal antibody (Native) (STJA0040583)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/ICC/WB
Host / Source	Chicken

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	10 mg/mL
Conjugation	Unconjugated
Purification	Mixture of IgY fraction and affinity-purified antibodies
Dilution Range	WB 1:5000-1:10000 WB Brain 1:5000-1:10000 IHC 1:2500-1:5000 ICC 1:2500-1:5000
Formulation	Sodium phosphate (10 mM, pH 7.2) buffered isotonic saline (0.9%, w/v) , glycerol (50%, v/v) , with thimerosal (0.01%, w/v) as an anti-microbial agent.
Isotype	IgY
Molecular Weight	117 kDa
Storage Instruction	Store at -20°C in the dark. Under these conditions, the antibodies should have a shelf life of at least twelve months, provided they remain sterile.

TARGET INFORMATION

Immunogen Region	Native
Immunogen Sequence	E Coli

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



IHC staining rat basolateral amygdala LacZ labeling using antibody green STJA0040583