

## Anti-Doublecortin antibody (STJA0040589)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC/ICC
Host / Source	Chicken
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	200 µg/mL
Conjugation	Unconjugated
Purification	Affinity-Purified IgY
Dilution Range	WB Not Recommended WB Brain Not Recommended IHC 1:1000-1:2000 ICC 1:1000-1:2000
Formulation	Phosphate-buffered (10 mM) isotonic (0.9%, w/v) saline ("PBS, pH 7.2) with sodium azide (0.02%, w/v) added as a preservative."
Isotype	IgY
Molecular Weight	40 kDa
Storage Instruction	Store at 4°C in the dark. Under these conditions, the antibodies should have a shelf life of at least twelve months, provided they remain sterile. For longer term storage, aliquot and freeze to avoid freeze-thaw of the antibody.

### TARGET INFORMATION

Gene ID	<a href="#">9201</a>
Gene Symbol	<a href="#">DCLK1</a>
UniProt ID	<a href="#">DCLK1_HUMAN</a>
Immunogen Sequence	Human, Mouse, Rat

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------