

Anti-V5 antibody (STJA0040595)

GENERAL INFORMATION


Product Type	Primary antibodies
Applications	ICC/WB
Host / Source	Chicken

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.2 mg/mL
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	WB 1:500-1:1000 ICC 1:100-1:500
Formulation	Phosphate-buffered (10 mM) isotonic (0.9%, w/v) saline ("PBS, pH 7.2) with sodium azide (0.02%, w/v) added as a preservative."
Isotype	IgY
Storage Instruction	Store at 4°C in the dark. Under these conditions, the antibodies should have a shelf life of at least twelve months, provided they remain sterile. For longer term storage, aliquot and freeze to avoid freeze-thaw of the antibody.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



IHC staining mouse retinal horizontal cell using Anti-V5 antibody with DAPI counterstain STJA0040595