

Anti-Giantin antibody (STJA0040602)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC
Host / Source	Chicken
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	200 µg/mL
Conjugation	Unconjugated
Purification	Affinity-Purified IgY
Dilution Range	IHC 1:100-1:500 ICC 1:100-1:500
Formulation	Phosphate-buffered (10 mM) isotonic (0.9%, w/v) saline ("PBS, pH 7.2) with sodium azide (0.02%, w/v) added as a preservative."
Isotype	IgY
Molecular Weight	376 kDa
Storage Instruction	Store at 4°C in the dark. Under these conditions, the antibodies should have a shelf life of at least twelve months, provided they remain sterile. For longer term storage, aliquot and freeze to avoid freeze-thaw of the antibody.

TARGET INFORMATION

Immunogen Sequence	Human
Specificity	No cross-reactivity reported

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



anti-Giantin antibody ICC staining Golgi apparatus in HeLa cells with green fluorescence
STJA0040602